

Enjoy water activities with the Neuro 2 Swim Kit



Robust and fully waterproof

With the Swim Kit, Neuro 2 users can enjoy a wide range of water activities – from swimming and playing in the pool, to sailing and rafting. The fully waterproof and reusable silicone cover is worn over the Neuro 2 sound processor and lets users get completely wet, without any concern for their Neuro 2.

Enjoy the water

Tailor-made for Neuro 2, the Swim Kit is useful for users of all ages in many situations. It enables children taking swimming lessons to hear the instructor, outdoor enthusiasts to take part in water activities and parents to be able to communicate with their child in the water. Wherever it is used, the Swim Kit allows users to still feel included and connected to their surroundings during water activities.

For saltwater, freshwater, chlorinated pools or soapy bathwater

The Swim Kit is extensively tested to stand up to many different situations, like all devices and accessories from Oticon Medical. The robust solution can be used in all types of water, for example, in the sea, in swimming pools, in freshwater rivers, or at home in the bath or shower. It also provides great protection against dust and humid environments.

The Swim Kit provides full waterproof protection at depths of three metres underwater for up to two hours and is designed to withstand exposure to:

- Sunlight
- Sweat
- High humidity levels
- Sun lotions
- Shampoos and conditioners



The difference is in the detail

The highest level of craftsmanship has gone into each and every detail of the Neuro 2 sound processor. It is designed to be thin and small, but at the same time robust and reliable. The Neuro 2 sound processor meets the IP68 standard for water and dust protection – the highest classification against water and dust ingress.

This high level of protection ensures that the sound processor will not be damaged if it is accidentally exposed to water at depths of up to one metre for 30 minutes.



Get wet with confidence

To keep the Swim Kit in place – even during lively activities in water – there are three different retention solutions included with each Swim Kit. These are designed to help prevent the Swim Kit from being dropped or lost.

Ear Plug Retention

The Ear Plug Retention solution is placed inside the ear for extra security. This instant-fit solution is available in four standard sizes for both right and left ears.

Retention Tubing

The retention tubing comes in clear material for discretion and is adjustable for all ear sizes. It can be used with all types of bathing suits, as it does not have to be attached to any clothes.

Thin Safety Line

This safety line is a good solution for one-piece bathing suits and rash vests and is made in a transparent material for discretion.



About the Swim Kit

- Compatible with both small and large rechargeable battery types
- One set includes two swim sleeves
- Each swim sleeve can be used for one year or 50 times, whichever comes first
- Comes with a practical compact storage case
- Available in 10 different variants for different combinations of magnet strength and antenna type
- Three standard retention solutions included

Tell your clients

- Enjoy the water with the fully waterproof Neuro 2 Swim Kit
- The Swim Kit offers full protection up to three metres underwater for up to two hours
- This robust solution can withstand exposure to sunlight, sweat, high humidity levels, sunscreen, shampoo and conditioner – and all types of water
- It can be used by all ages



Visit [oticonmedical.com](https://www.oticonmedical.com) to see the video on how to use the Swim Kit